

MiDSystem Tutorial

Table of Contents

| | |
|---|----|
| General Instruction..... | 2 |
| Task Submission | 3 |
| Data upload | 3 |
| Parameter settings | 4 |
| Phylogenetic Tree | 9 |
| Task Status..... | 9 |
| Result Pages..... | 10 |
| Overview..... | 10 |
| Quality Control, De Novo Assembly, Gene Prediction, and GO Term Annotation..... | 11 |
| Phylogenetic Tree | 13 |
| Result pages of the metagenomics pipeline | 14 |
| System Info | 14 |

General Instruction

MiDSys^{tem} consists of four pipelines: reference-guided genome assembly, non-reference-guided assembly, long-read/hybrid assembly, and metagenomics analysis.

Reference-guided genome assembly and non-reference-guided assembly pipelines are for the single microbial species with short-read sequencing data (e.g., Illumina). If you already have a reference genome for your assembled sequences, we suggest conducting the reference-guided pipeline. In addition, MiDSys^{tem} provides assembly pipelines of the single microbial species with long-read sequences (e.g., Pacific Biosciences and Oxford Nanopore Technologies). Also, the hybrid assembly allows performing genome assembly with short-read and long-read sequences simultaneously. Before submission, please make sure what kind of experimental design in the uploaded files and select a suitable pipeline for analysis. **Trying to run metagenomics samples on these pipelines above, or to run single species data on the metagenomics pipeline will fail the analysis process** since they are designed for different purposes. The pipelines menu locates at the top of the home page (Fig. 1).




Fig. 1. The menu at the top of the home page allows users to select the pipelines to run.

There are step-by-step instructions on each pipeline's home page. Please follow them. For all the pipelines, you have to input your basic information first (i.e., the email address) (Fig. 2), so that MiDSys^{tem} can send a unique link to check the status or the final report of the submitted task. The unique link for your task is generated for security reasons. With the correct email address and link, you can come back to view your data at any time, so feel free to close the page after submission. Notification mails will also be sent to the email address when the analysis process is finished (either success or failure).

The image shows a web form titled 'Basic Information' in a light red header. Below the header, there is a text prompt 'please provide your basic information.' followed by a label 'Email address:' and a text input field. Below this, there is another text prompt 'Please re-enter your Email for confirmation.' followed by a label 'Re-enter Email address:' and a second text input field.


Fig.2. MiDSys^{tem} will send a unique link for your task according to the given email address.

Task Submission

Data upload

Data upload is required for all the pipelines. The R1/R2 and/or long-read files should be uploaded separately in **.fastq/.fq**, and **.gz** formats. The maximum size for each file is **20 GB**. Users can choose to upload through their browser (Fig. 3A) or with URLs (e.g. google drive links)(Fig. 3B). The reference-guided pipeline will require an extra reference file to upload in **.fasta/.fa**, **.fna**, and **.gz** formats.

(A)



Step 1. Reads Upload

Select a method for uploading files:
☒ From browser ☐ From URLs

In this step, we will provide FastQC, MultiQC and Trimmomatic to control the quality of your data uploaded below.

Please upload your sequencing files here. Only allow **.fastq/.fq** and **.gz** format. Allowed maximum size for each file: **20 GB**.

[+ Add files...](#)

1). demo_R1.fastq.gz (540.6 MB)

☒ R1 ☐ R2

[Cancel](#)


2). demo_R2.fastq.gz (540.7 MB)

☐ R1 ☒ R2

[Cancel](#)

[Start upload](#)

(B)



Step 2. Reference Upload

Upload a bacterial reference genome for the reference guided assembly.

Please provide a URL of your reference file here. Using a share link from Google Drive is available. Only allow **.fasta/.fa**, **.fna**, and **.gz** format.

Reference:

[Confirm the URL](#)

Fig. 3. MiDSytem provides two upload methods for data submission. (A) Upload data from browser. (B) Upload data from an existing URL.

In Long-Read/Hybrid assembly, you can consider further uploading short-read sequences after uploading long-read sequences if available (Fig. 4). The long-read and hybrid assembly pipelines are automatically switched depending on the short-read files are provided or not.

Step 2. Short-Read Upload

Does upload short-read sequences for the hybrid assembly?

☐ No
☒ Yes, use the hybrid assembly.

In this step, we will provide FastQC, MultiQC, and Trimmomatic for short-read sequences to control the quality of your data uploaded below.

Please provide URLs (e.g., <http://>, <https://>, or <ftp://>) of your sequencing files in R1 and R2 fields. Using a share link from Google Drive is available. Only allow .fastq/.fq and .gz format.

R1:

R2:

Fig. 4. Upload short-read files for conducting the hybrid assembly pipeline in Long-Read/Hybrid assembly.

Parameter settings

After file upload step, the following analysis settings are enabled. Default settings are provided to all the tools that we will run in the analysis process (Fig 5). Users can also customize the parameters based on their needs. We provide default settings, but you can use the “Customized” option to adjust the parameters (Fig 6).

Step 2. De Novo Assembly

In this step, we provide A5-miseq and several tools for assessment of the assebled sequence.

☒ Default Settings
☐ Customized

Step 3. Gene Prediction

In this step, please select one of the gene prediction tools below:

☒ GeneMark
☐ Augustus

Step 3(conti.). Predicted Gene Assessment

In this step, we provide BUSCO for assessment.

☒ Default Settings
☐ Customized

Step 4. GO Term Annotation

In this step, we provide InterProScan for GO term annotation

☒ Default Settings
☐ Customized

Fig. 5. Default settings of non-reference assembly.

Step 2. *De Novo* Assembly

In this step, we provide A5-miseq and several tools for assessment of the assebled sequence.

☐ Default Settings
☒ Customized

| Quast settings | Values |
|---------------------------------|--------|
| minimum contig-thresholds (>=0) | 300 |

| BUSCO settings | Values |
|----------------|------------------|
| species | Escherichia coli |
| e-value | 1e-03 |

| Bowtie2 settings | Values |
|--|---|
| --no-unal (Suppress SAM records for reads that failed to align) | <input type="radio"/> No <input checked="" type="radio"/> Yes |

Fig 6. Customized parameters for relevant tools in the *de novo* assembly step.

The parameter settings would be different depending on the selected pipeline, but the parameters that you can adjust are displayed after you choose “Customized” option. Here, we are listing parameters of assembly pipelines in MiDSystem that can be customized and their descriptions (Table 1).

Table 1. Parameters of assembly pipelines in MiDSysstem.

| Parameter | Affected Tool | Assembly Pipeline | Default Value | Data Type | Description |
|---------------------------|---------------|--|------------------|--|--|
| mode | Unicycler | Long-read, Hybrid | Normal | {Conservative, Normal, Bold} | Conservative mode is least likely to produce a complete assembly and smaller contigs but has a very low risk of misassembly. Bold mode is most likely to produce a complete assembly and longest contigs but carries greater risk of misassembly. Normal mode is intermediate regarding both completeness and misassembly risk, as well as moderate contig size. |
| minimum contig-thresholds | Quast | Reference-guided (for contigs), Non-reference, | 300 | Non-negative integers | Lower threshold for contig length. |
| | | Reference-guided (for scaffolds), Long-read, Hybrid | 1000 | | |
| species | BUSCO | Reference-guided, Non-reference, | Escherichia coli | {Escherichia coli, Staphylococcus aureus, | Name of existing Augustus species gene finding parameters. |

| | | | | | |
|-------------|--------------|---|------------------|---|---|
| | | Long-read, Hybrid | | thermoanaerobacter_tengcon gensis} | |
| e-value | BUSCO | Reference-guided, Non-reference, Long-read, Hybrid | 1e-03 | Positive floats | E-value cutoff for BLAST searches. |
| --no-unal | Bowtie2 | Reference-guided, Non-reference, Hybrid | Yes | {Yes, No} | Suppress SAM records for reads that failed to align. |
| species | Augustus | Reference-guided, Non-reference, Long-read, Hybrid | Escherichia coli | {Escherichia coli, Staphylococcus aureus, thermoanaerobacter_tengcon gensis} | Select a training set of specific species to predict genes. |
| strand | Augustus | Reference-guided, Non-reference, Long-read, Hybrid | Both | {Both, Forward, Backward} | Predict genes on which strands. |
| e-value | BLAST | Reference-guided, Non-reference, Long-read, Hybrid | 1E-5 | Positive floats | The number of hits one can “expect” to see by chance when searching a database of a particular size. |
| file format | InterProScan | Reference-guided, Non-reference, | tsv | {tsv, xml, gff3, svg, json, html} | Select output formats of InterProScan reports. |

| | | | | | |
|----------------------|----------|----------------------|------|-----------------------------|---|
| | | Long-read, Hybrid | | | |
| output file format | GraPhlAn | Metagenomics | png | {png, pdf, ps, eps, svg} | Select output formats of cladograms. |
| clustering threshold | CD-HIT | Metagenomics | 0.97 | Positive floats (0.7-1.0) | The of threshold sequence identity. This is the default cd-hit's "global sequence identity" calculated as: number of identical amino acids in alignment divided by the full length of the shorter sequence. |
| insert size (-X) | Bowtie | Metagenomics | 800 | Positive integers (250-800) | Maximum insert size for paired-end alignment. |
| mismatch (-v) | Bowtie | Metagenomics | 3 | Non-negative integers | Report end-to-end hits with <= this number of mismatches. |
| e-value | HMMER | Metagenomics | 1E-5 | Positive floats | Report sequences <= this e-value threshold in output. |

Phylogenetic Tree

The phylogenetic tree step is available in all assembly pipelines of single species (Fig. 7). It is optional, select “Yes” for construction of a phylogeny with your data and other 10 species. The database will be loaded in to the page and the selections of species will be enabled. (Note: If the loading is failed, try to choose “No”, then “Yes” again to reload the data.) You can name your sample with at most 10 characters, only A-Za-z0-9 and _ are allowed. The default value is “my_sample” (i.e., the tree will display my_sample to represent your sample).

Fig. 7. Setting for construction of a phylogenetic tree.

The “Search” box helps find out the species of interest. The species that you selected will be displayed on the right panel (Fig. 8). (Note: when you select more than 10 species, the oldest ones that you selected before will be removed.)

Fig. 8. Selection species in the MiDSysSystem database for construction of the phylogeny with uploaded data.

Task Status

After submission, you will receive an email with a link to help check out the status of the submitted task Fig. 9. A certain step will display “SKIPPED” when they are not necessary for the specific pipeline. After the task is completed, another email will be sent to the mail address

you provided, and this page will be no longer available since it will be automatically redirect to a report page.




Fig. 9. Screenshot of the task status.

Result Pages

Overview

The result page is shown in Fig. 10. The name of performed pipeline displays at the top left corner. The task information and download link of full analysis reports are provided for users. The left panel can be used to navigate through summary reports of tools. To make complicated reports that can be interpreted, all outputs in MiDSysyem are summarized as tables, pie charts, bar charts, etc.




Fig. 10. Screenshot of the result page.

Quality Control, De Novo Assembly, Gene Prediction, and GO Term Annotation

Most of the visualizations are intuitive. To sum up, the “Quality Control” page includes detailed information about the submitted and trimmed R1/R2 sequences. The “De Nove Assembly” and “Reference-Guided Assembly” pages provide N50 and the quality assessment reports for assembly processes. Gene prediction results are also provided with visualized assessment reports. For functional annotation, the GO term annotation bar chart is displayed corresponding to three categories.

The dashboard on the left panel will be modified based on the pipeline selected. For example, the Reference-Guided Assembly Result provides summary statistics about the assembled genome after construction of longer scaffolds from short contigs (Fig. 11), but the page link does not show on the dashboard of other pipelines. Some original reports analyzed by specific tools can be directly downloaded via provided links.




Fig. 11. The Reference-Guided Assembly Result page shows the summary statistics of scaffolds produced by Ragout. The assembled genome can be downloaded by clicking the download button.

Three categories of the GO term annotation are displayed as bar charts in different colors (Fig. 12), i.e., biological process (green), cellular component (orange), and molecular functions (blue). The GO term bar charts provide for the *de novo* assembly pipelines of single species are based on frequency plots (i.e., how many times the specific GO term exists). However, in the metagenomics pipeline, the GO term bar charts are based on the z score. Only the top-ranked GO terms, either based on frequency or z score, for each category are displayed on website pages. Detail about z score calculation is provided in the Method section of the manuscript.




Fig. 13. The bar charts show z scores of significant GO terms in the metagenomics pipeline.

Phylogenetic Tree

Phylogenetic tree construction is only provided for the single species pipelines. A result of the *E. coli* EC4437 strain sample is provided in the manuscript case study. NCBI external links for each species selected are available on the result page (Fig. 14).




Fig. 14. Screenshot of the phylogenetic result.

Result pages of the metagenomics pipeline

On the taxonomic abundance page of the metagenomics results, a cladogram, bar chart, and table are used to present the species composition in a sample. The protein domain information page is available as Fig. 15. The domain frequency and description are displayed in Frequency Table. The frequency column is the number of genes mapped to the specific domain. Accession numbers are external links to the Pfam database. This table can be downloaded in a CSV file.

Metagenomics Pipeline

Dashboard

Overview

Quality Control

Assembly

Taxonomic Abundance

Gene Prediction/Clustering

Functional Annotation and Abundance

Domain Mapping

System Info

Frequency Table

Show 10 entries

Search:

| Domain | Accession(Link) | Frequency | Description |
|--------------|----------------------------|-----------|---|
| 2-Hacid_dh_C | PF02826.18 | 2 | D-isomer specific 2-hydroxyacid dehydrogenase, NAD binding domain |
| 4HBT | PF03061.21 | 1 | Thioesterase superfamily |
| 5_3_exonuc_N | PF02739.15 | 1 | 5'-3' exonuclease, N-terminal resolvase-like domain |
| AAA_16 | PF13191.5 | 2 | AAA ATPase domain |
| AAA_18 | PF13238.5 | 1 | AAA domain |
| AAA_21 | PF13304.5 | 2 | AAA domain, putative AbiEii toxin, Type IV TA system |
| AAA_22 | PF13401.5 | 1 | AAA domain |
| AAA_31 | PF13614.5 | 1 | AAA domain |
| AAA_33 | PF13671.5 | 1 | AAA domain |
| ABC-3 | PF00950.16 | 1 | ABC 3 transport family |

CSV

Showing 1 to 10 of 617 entries

Previous12345...62Next

Fig. 15. Frequency Table shows on the domain mapping page.

System Info

The link of System Info can be found on the dashboard in all pipelines. This page shows information about the MiDSystem version and collected all third-party tools and databases with their versions and executing functions that were conducted in the current pipeline (Fig. 16). This information would help the reproducibility of analysis results from MiDSystem. The names of tools and databases on the tables contain external links to their official websites.

Long-Read Assembly

Dashboard

Overview

Quality Control

De Novo Assembly

Gene Prediction

GO Term Annotation

Tree

System Info

System Info

The page provides information about the MIDSystem version, lists of third-party tools, and databases that were executed in this task.

MIDSystem Version

20210610

List of Third-party Tools

| Tool Name | Version | Executing Function |
|------------------------------|-----------|--|
| BLAST+ | 2.6.0 | Genome Assessment, Gene Annotation, Gene Prediction Assessment |
| BUSCO | 2.0 | Genome Assessment, Gene Prediction Assessment |
| ETE 3 | 3.1.1 | Phylogenetic Tree |
| GeneMarkS | 4.32 | Gene Prediction |
| InterProScan | 5.25-64.0 | Functional Annotation |
| LongQC | 1.2.0b | Quality Control |
| MUSCLE | 3.8.1551 | Phylogenetic Tree |
| OrthoMCL | 2.0.9 | Phylogenetic Tree |
| QUAST | 4.4 | Genome Assessment |
| RAxML | 8.2.10 | Phylogenetic Tree |
| Racon | 1.4.21 | Genome Assembly |
| Unicycler | 0.4.9b | Genome Assembly |

List of Databases

| Database Name | Released Version/Date |
|---|-----------------------|
| BUSCO bacteria odb9 | v2 |
| Ensembl bacteria database | Release 32 |
| InterPro | v64.0 |
| RefSeq non-redundant proteins | 03/13/2018 |

Fig. 16. Screenshot of the System Info page.